



RAPID ALCOHOL DISINFECTION

Incidin™ Liquid

ALCOHOL BASED READY-TO-USE SURFACE DISINFECTION

In addition, it offers high reliability to the user through its broad spectrum of efficacy including TB and Noro.

Incidin™ Liquid also provides:

- Short contact times
- Dye and fragrance free formulation
- Excellent material compatibility on alcohol resistant surfaces
- Compatibility with Incidin™ Dry Wipes System



Incidin™ Liquid

Alcohol based ready-to-use surface disinfection



Incidin™ Liquid
650 ml Bottle

Incidin™ Liquid
1 L Bottle

Incidin™ Liquid
5 L Canister

BROAD EFFICACY SPECTRUM IN SHORT CONTACT TIMES

Due to its broad efficacy spectrum in short contact times, Incidin™ Liquid is ideal even for intensive care area and other high risk areas. Contact times are detailed in table below.

RELIABLE USE

Being ready-to-use and offering proven material compatibility on alcohol-resistant surfaces, Incidin™ Liquid ensures the reliable use for all application areas.

DYE AND FRAGRANCE FREE

The dye and fragrance free formulation leaves no active residues and ensures a convenient product use.

MICROBIOLOGICAL AND VIROLOGICAL TEST RESULTS

EFFICACY UNDILUTED	CONTACT TIME
Disinfection of surfaces acc. to VAH with / without mechanical action (bactericidal, yeasticidal, tuberculocidal)	5 min
Rapid surface disinfection acc. to VAH method 2015 without mechanical action	1 min
EN 13727 / EN 13624 / EN 14348 (bactericidal, yeasticidal, tuberculocidal)	1 min
EN 16615 bactericidal / yeasticidal	1 min
EN 16615 (method) tuberculocidal	1 min
RKI-DVV 01/2004* / EN 14476 enveloped viruses	30 sec
DVV 2015* / EN 14476 limited spectrum viricidal	10 min
EN 14476 Adeno	1 min
EN 14476 Noro	10 min
DVV 2015* / EN 14476 (clean) Rota	30 sec
EN 16777 Noro (dirty)	10 min
EN 16777 Adeno (dirty)	5 min
EN 16777 Vaccinia (dirty)	1 min

*DVV 2015 (suspension test)

ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Incidin™ Liquid Bottle	12 x 650 ml	XXXXXX
Incidin™ Liquid Bottle	12 x 1 L	XXXXXX
Incidin™ Liquid Canister	2 x 5 L	XXXXXX

APPLICATION AREA

Fast acting, ready-to use surface disinfection solution for medical devices.

INSTRUCTIONS FOR USE

Use product undiluted on a non-woven wipe. Wipe the alcohol resistant surface and the product work, or spray surface from ca. 30cm distance (max 40ml/m²) and let dry. Ensure the entire surface is completely wet.

INGREDIENTS

100g contains the following:
35g 2-Propanol, 25g 1-Propanol

STORAGE

Product should be stored in the original container below 25°C in a cool, dry location out of direct sunlight.

MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.

ECOLAB EUROPE
Richtstrasse 7
8304 Wallisellen
Switzerland



ECOLAB DEUTSCHLAND GmbH
Ecolab-Allee 1
40789 Monheim am Rhein
Germany

ECOLAB®

www.ecolab.com