

# Klerwipe™ Sealed Edge Knitted Dry Wipes

PHARMACEUTICAL

General purpose cleanroom wipe. Ultrasonically sealed edge wipe is 100% knitted polyester and laundered and packaged in an ISO Class 4 cleanroom. Recommended for use in ISO Class 3 and above environments.

The validated sterile version of this wipe is excellent for cleaning applications during aseptic processing activity including Engineering Controls (PEC) such as laminar air hoods, biological safety cabinets and compounding aseptic isolators where sterility of end product must be assured.

## Knitted

### FEATURES

- Soft hand and edge for sensitive surfaces, abrasion resistant
- Low particle and fiber generation, high level of cleanliness
- Chemically compatible with disinfectants and cleaning solutions
- Sterile version available
- Sealed edge, cleanroom laundered, knitted wipes offer higher purity and durability compared to nonwoven wipes

## Validated Sterile

### FEATURES

- Gamma irradiated to a 10<sup>-6</sup> Sterility Assurance Level (SAL)
- Certificates of Sterility available
  - Expiration date
  - Radiation dosage
  - Lot information
  - Meets cGMP requirements for traceability
- Double bagged
- Header style packaging
- Package Integrity testing



## ORDER INFORMATION

### KNITTED

PART #	53000349
SIZE	9X9"
WIPES/PACK	150
PACKS/CASE	20

PART #	53000352
SIZE	12X12"
WIPES/PACK	150
PACKS/CASE	10

### VALIDATED STERILE

PART #	53000347
SIZE	9X9"
WIPES/PACK	10
PACKS/CASE	60

PART #	53000348
SIZE	9X9"
WIPES/PACK	20
PACKS/CASE	40

PART #	53000350
SIZE	12X12"
WIPES/PACK	10
PACKS/CASE	40

PART #	53000351
SIZE	12X12"
WIPES/PACK	20
PACKS/CASE	25

FOR ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

### WORLDWIDE HEADQUARTERS

1 Ecolab Place  
St. Paul, MN 55102  
USA

[www.ecolab.com/lifesciences](http://www.ecolab.com/lifesciences)  
[www.ecolab.com/pharmaceutical](http://www.ecolab.com/pharmaceutical)



Everywhere It Matters.™