



AseptoClean

Chlorinated alkaline cleaner for milking equipment and milk tanks.

- ▲ Good performance in cleaning
- ▲ Cost-efficient solution

APPLICATION

Pre-clean with water. Clean with **AseptoClean** at 0.5% (minimum 10 min. at 40°C - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithClean**.



FEATURES AND BENEFITS

AseptoClean is an alkaline chlorinated product with boosted cleaning properties. It is suitable for all kinds of standard milking equipment, including cooling tanks.

The **Clean** products are the most applicable basic program for good milking machine hygiene.

SPECIFICATIONS

AseptoClean	
Physical Form	Light yellow liquid
Active Substance	Sodium hydroxide, Sodium hypochlorite
pH at 20°C:	13 - 14
Density at 20°C (g/cm ³):	1.19 - 1.22
Storage temperature:	-5°C to +30°C

LOGISTICS

Packaging	12 kg	24 kg	72 kg	245 kg	1175 kg
Product code	2321510	2321530	2334630	2321550	2321560
EAN code	4028162321517	4028162321531	4028162334630	4028162321555	4028162321562

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com