


HorolithClean



Acidic cleaner for milking equipment and milk tanks.

- ▲ Good performance in descaling
- ▲ Balance of two acids
- ▲ Cost-efficient solution

APPLICATION

Pre-clean with water. Clean with **HorolithClean** at 0.5% (minimum 10 min. at 40°C - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **AseptoClean**.



FEATURES AND BENEFITS

HorolithClean is an acidic product which assures very good scale removing and cleaning properties. It is suitable for all kind of standard milking equipment, including cooling tanks. The **Clean** products are the most applicable basic program for good milking machine hygiene.

SPECIFICATIONS

HorolithClean	
Physical Form	Clear, colourless liquid
Active Substance	Inorganic acids
pH at 20°C:	1.0
Density at 20°C (g/cm ³):	1.20 - 1.22
Storage temperature:	-30°C to +40°C

LOGISTICS

Packaging	12 kg	24 kg	72 kg	250 kg	1200 kg
Product code	2331990	2332000	2334740	2332020	2332030
EAN code	4028162331998	4028162332001	4028162334746	4028162332025	4028162332032

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com