



# Des Sanitary

## WASHROOM CLEANER DISINFECTANT

POWERED BY  
**ECOLAB**



### DESCRIPTION:

A daily washroom cleaner disinfectant with excellent cleaning performance to remove body fats, soap residues and dirt on hard surfaces.

- Provides excellent detergent action whilst effectively removing limescale
- Rapid broad spectrum antimicrobial performance for improved safety in the washroom environment
- Maximum cleaning performance leaving a pleasant fragrance



### USE INSTRUCTIONS:

#### FOR HARD SURFACES AND FLOORS

1. Always wear PPE when handling the product, place warning signs out for safety
2. Dilute 5 - 30 ml of Sanitary per 1L of water, spray bottle dilution at 20-30ml
3. Apply solution to surface dispensing via spray bottle or bucket
4. Wipe and clean surfaces and leave in contact for the required time shown
5. Rinse with water or wipe with a clean damp cloth and allow to air dry

#### FOR TOILET BOWL CLEANING

1. Always wear PPE when handling the product, place warning signs out for safety
2. Use product neat via 1L bottle
3. Close toilet lid and flush toilet, open toilet lid and reduce water level by pushing water repeatedly with brush
4. Apply product under rim, around toilet bowl and under water line
5. Use toilet brush to ensure all surfaces are covered and to remove dirt, where possible leave brush in toilet
6. Leave in contact for 5 minutes maximum, flush and rinse toilet and brush



### EFFICACY:

Bactericidal acc. to EN 1276 dirty conditions at 1.5% for 5 minutes
Bactericidal acc. to EN 1276 clean conditions at 1% for 1 minute
Bactericidal acc. to EN 16615 dirty conditions at 2.5% for 1 minute
Yeasticidal acc. to EN 16615 dirty conditions at 2% for 1 minute
Virucidal acc. to EN 14476 enveloped viruses dirty conditions at 2% for 5 minutes
Virucidal acc. to EN 14476 enveloped viruses clean conditions at 3% for 1 minute
HEALTHCARE USE ONLY
Yeasticidal acc. to EN 1650 dirty conditions at 2% for 5 minutes
Yeasticidal acc. to EN 13697 dirty conditions at 7% for 15 minutes



### APPLICATION:

Our certified disinfectant products are 99.99% effective against bacteria and enveloped viruses found in public spaces when used with the correct cleaning procedure.

For all acid-resistant surfaces and floors in the washroom, not to be used with products containing chlorine.

For easy toilet and urinal cleaning the 1L bottle is fitted with a directional nozzle.



### SAFETY & STORAGE:

- Do not mix with any other product. Use cleaner disinfectant safely. Always read the label and product information before use. Safety data sheet available for professional users on request.
- Dispose of packaging only if emptied and closed. Product residues: See Safety Data Sheet.
- Active substance 100g liquid contains 8g L-(+)-lactic acid, 2,45g Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-C16)).
- Store between -5 °C and 40 °C in original closed container only.

PRODUCT CODE:    PACK SIZE:

9104740            2x5L

9106550            12x1L



### ADDITIONAL RESOURCES



- MSDS
- WALLCHARTS
- PRODUCT GUIDES



Worldwide Headquarter  
Ecolab Inc.  
1 Ecolab Place, St. Paul, MN55102  
1.800.35.CLEAN

European Headquarter  
Ecolab Europe GmbH  
Richtistr. 7  
CH-8304 Wallisellen